

FOFO OVER **20** YEARS EXPERIENCE

IN HIGH-FREQUENCY PRODUCTS



COMPANY INTRODUCTION

Foshan Hongfeng Co., Ltd mainly produces medical rehabilitation products, specializing in the design and manufacture of air compression therapy machines, compressor nebulizers, suction machines, air cushions and other rehabilitation products. Foshan HongFeng has over 20 years of experience in making high-frequency products. Our high frequency machines are 225KW, and are the most technologically advanced with professional patented designs. There are over 360 workers in the group. Our company has nstrong R & D capabilities, from mold design and manufacturing to software research and development, high-frequency product design and manufacturing are all completed internally, and quality control meets European and American requirements.

Our products have been approved by FDA, 510K, ISO and CE MDR. The sales team's languages include English, Russian, Spanish, French, and Japanese. We have exported to the USA, Europe, Asia, the Middle East, and many other countries in the world. FOFO has earned a great reputation in our industry amongst our customers and peers. We sincerely hope we can build a stable business relationship with you and welcome you to visit our factory. We guarantee high-quality products, competitive prices, prompt delivery schedules and reliable services.



CE 0598 FDA ISO13485 510K



FOFO



FOFO

Air Compression Therapy Systems



Foshan Hongfeng Co.,Ltd

Address: No.4-2 Leqiang Road, Leping Sanshui,Foshan,528100,GuangDong,China
Tel: 0757-88791710 ext 888 0592-2097262
Email: fofo@fomedical.com
Mob: +86 13055506651
www.fomedical.com



YoYube



www.fomedical.com

Air Compression Therapy System Cuffs




APPLICATION

Lymphedema treatment typically includes a few methods and compression devices is one of them. The pneumatic compression pump provides on-and-off pressure to keep fluid moving through lymph vessels and veins so it doesn't build up in arms, legs and other parts of your body. The device works by connecting the pump to a sleeve that wraps around affected areas. The pump drives a regular cycle of pressure that inflates and deflates the sleeve. Compression pumps may also reduce complications of lymphedema.

How can air compression therapy systems help ?

- ✓ Improving blood and lymph circulation.
- ✓ Pushing the residual lymph towards the lymph nodes.
- ✓ Removal of lymphatic edema and exudates.
- ✓ Improving the nutrition and oxygenation of the skin and other tissues.




Sleeve	Description	Size
Leg	Leg Cuff M	91*65cm
	Leg Cuff L	100*70cm
	Leg Cuff XL	110*70cm Overlapping
	Leg Cuff XXXL	125*76cm Overlapping
Arm	Arm Cuff M	87*51cm
	Arm Cuff L	90*60cm Overlapping
Waist	Waist Cuff M	120*37cm
	Waist Cuff L	130*37cm Overlapping
Trousers	Short Pants (4,6,8 Chambers)	80*76cm
	Long Pants (6,8 Chambers)	142*85cm
Vest/cloth	Vest (6,8 Chambers)	78*72cm
Jacket	Double Arm/Shoulder	133*37 (pair)cm
hand	Left/Right	M/L/XL



ABC DEF 6 Modes



- A mode** One by one, single part massage, after all chambers are treated, one cycle is complete
- B mode** focus therapy from bottom to top
- C mode** Therapy is focused on massage from bottom to top
- D mode** All chambers inflate and deflate together
- E mode** Combined Sequential Mode
- F mode** Combined Focused Mode







Wireless Magnetic Portable Air Compression Boots for Travel

- System consists of a detachable electronic main unit and garments.
- The controller is attached to the magnetic panel and can be removed.
- Allows multi-use of attaching the main unit to legs, arms and waists,pants.

Model	FO3003D
Device Size	12.5x7x4.8cm
Massage Modes	4 modes
Pressure Range	80 ~ 250 mmHg
Treatment Time	10 ~ 60 mins
Chambers	5
Input	DC5V/3A/2600mA
Cuffs Compatibility	Leg, Arm, Waist, Long Pants



LCD touch screen Air Compression Therapy System with Bluetooth



SPECIFICATIONS

- Multi-Lingual Device
- Demo Videos
- Power Consumption: 30W/40VA
- Device Size: 25cm*12cm*9cm
- Massage Mode: 6 modes
- Pressure Range: 30~270 mmHg
- Treatment Time: 10~60 mins
- Input: DC 15V/5200mA
- Each Chamber can be Partially Skipped
- Can use as USB Mobile Phone Charger
- Weight: 1.4 kg



FO3002E



App Optional



7-inch large screen



APP synchronization control



Pre set cold therapy mode



2.5-4 hours of battery life



USB battery external discharge



Habitual memory

WiFi and Bluetooth Connectable/ App Control


Built-in Video Instructions

Full-Body Garments


Intermittent Pressure Compression Device (IPC)

4/6/8 Chambers Medical CE FDA 510K

- Powerful Medical-Grade Air Pump Provides 250mmHg Pressure
- Upgraded A-F 6 Modes
- Chambers Optional Function skips sections which you don't want to compress
- Sleeves with overlapping design provide more comfort
- Extension Straps for wider Sleeves
- Easy-Split Designs for 2-4 Sleeves parallel use



Remote control distance <5 meters



FO3001-2



Specifications:

Pump	Sleeves	
Item: FO3001, FO3001-2	Leg(4/6/8 Chambers)	M 91*65cm / L 100*70cm
Power Supply: 110V/220V/50-60Hz		XL 110*70cm overlapping / XXXL 125*76cm overlapping
Air Output: 15-20 liter/min.	Arm(4/6/8 Chambers)	M 87*51cm / L 90*60cm overlapping
Pressure Range:30mmHg - 250 mmHg	Waist(4/6/8 Chambers)	M 120*37cm / L 130*37cm overlapping
Time: 10-60 min.	Long Pants (6/8 Chambers)	142*85cm / 142*76cm overlapping
Power Consumption: 45VA	Short Pants (4/6/8 Chambers)	82*77cm
Size: 24cm(L) x 20cm(W) x 11cm(H)	Jacket (6/8 Chambers)	102*73cm
Weight: 2.45kg	Vest(6/8 Chambers)	78*72cm

Sequential Compression Device (SCD)


4/6/8 Chambers MDR CE 510K

- Medical CE, FDA510K ClassII
- Li-ion Battery 2600mA
- Powerful Medical-Grade Air Pump Provides 250mmHg Pressure
- Upgraded A-E 5 Modes Instead of the old Two-Modes Massage and Therapy

FO3002



4.3" Screen Air Compression Therapy With Li-ion Battery



FO3002-1

Specifications:

Pump	Sleeves	
Item: FO3002 / FO3002-1	Leg(4/6/8 Chambers)	M 91*65cm / L 100*70cm
Power Supply: 100-220V/12V/2600mAH		XL 110*70cm overlapping / XXXL 125*76cm overlapping
Air Output: 15-20 liter/min.	Arm(4/6/8 Chambers)	M 87*51cm / L 90*60cm overlapping
Pressure Range:30mmHg - 250 mmHg	Waist(4/6/8 Chambers)	M 120*37cm / L 130*37cm overlapping
Time: 10-60 min.	Long Pants (6/8 Chambers)	142*85cm / 142*76cm overlapping
Power Consumption: 65VA	Short Pants (4/6/8 Chambers)	82*77cm
Size: 25cm(L) x 12cm(W) x 9cm(H)	Jacket (6/8 Chambers)	102*73cm
Weight: 1-1.6 kg	Vest(6/8 Chambers)	78*72cm

NEW

Premium 12 Chambers Air Compression Therapy System

Segmented Pressure Control:

- Each chamber can have the pressure adjusted, with sections from the ankle to the thigh.
- Creates a gradient pressure from the distal end to the proximal end.
- Promotes unidirectional flow of blood and lymph, and prevents stasis.

Dynamic Compression Mode:

- The airbags inflate and expand in sequence, creating a wavy compression.
- Similar to the effect of a "muscle pump".
- Accelerates the return of venous blood and alleviates swelling in the lower limbs.

F03012



Specifications:

Pump	
Item No.	F03012
Name	12-chambers lymphedema compression massage machine / compression therapy boots
Control Way	LED touch keypad and remote controller
Machine Size	30.2 × 25.4 × 15.8 cm
Pressure Range	30~250 mmHg
Treating Time	10-60 minutes (electronically adjustable in steps)
Programs	8 modes: A/B/C/D/E/F/G/H
Intensive/Skip Care	12 chambers can be individually switched on/off to target specific muscles, areas, or avoid injured/surgical sites


NEW

Hand-held 6 Chambered legs Compression Therapy

- F03003C is individual valves to the 6 chambers for legs with 5 modes. Max pressure 250mmHg perfect performance for compression.
- Portable handy 4 and 6 Chambers device easy for carry and storage.
- Small individual box is good for online shops.

Model	F03003C
Device Size	23.8*8.8*5cm
Weight	0.66KG
Chamber	4 and 6 chambersleeves
Modes	A/B/C/D/E
Pressure	30-250mmHg
Power Consumption	DC12V/2A
Maximum Power	30W
Noise Level	Below 55 dBA
Certificates	CE

F03003C



Heating+Compression therapy

CE 0598 FDA ISO13485


Pressurization and Heating Simultaneously Accelerate Blood Flow



F03001M-1/ F03002-2

SPECIFICATIONS

- Size: 23.5 x24x10.5cm
- Machine Voltage: AC 220V
- Power: 30W/40VA
- Weight: 2.2KG
- Heating Module voltage: DC 15V
- Feeling Temperature: approx. 40°C



High Quality Air Compressor



Heated Leg Straps



Targeted Therapy for Foot, Calf, Knee or Thigh



Effective for Lymphatic Edema



40°C Max

Cryo+Compression Therapy

Full Leg Cooling and Compression

SPECIFICATIONS

- Name: Cryo+Compression Therapy
- Dimension: 290*270*285mm
- Weight: 3.1KG
- Water Flow: 0.85L/MIN
- Volume: 6.5L
- Temperature: 5°C

F05001



Lower cuff Air Massager with heating

- F03003B is 3 chambers for legs. Max pressure 160mmHg .
- Small individual box is good for online shops.

SPECIFICATIONS

- Device Size: 20*6*4.5cm
- Weight: 0.65KG
- Chamber: 6 chamber sleeves
- Modes: F1/F2/F3/F4
- Pressure: 60-160mmHg
- Power Consumption: DC12V/1.5A
- Maximum Power: 14W
- Noise Level: Below 55 dBA
- Characteristics: Heating Pressure



F03003B



Portable DVT/Wireless Calf Massager

SPECIFICATIONS

- **Dvt:** Pressure 10-60mmhg , Sequential Compression, Intermittent Time 24/48/60 Sec, Type-c
- **Massager:** 3 Modes And Heating, Pressure 60/80/99mmhg, Heating Temperature 45°C-55°C

- Main Device Size: 14*7*0.25cm
- Input Voltage: DC 5V/2000 mA
- Carton Box: 380*340*380mm; 10pcs per carton; 5800pcs for 20ft, 13870 for 40HQ
- Certificates: ISO, FDA, Massage CE


F03003



Back stretching pad

Automatic Back Stretching Air Pad, Lumbar Spine Relaxer, Stretcher, Waist Tensioner, Loose Top, Easy-Roll up Mat, Back Exercises
Single pad 6/8 zones fit all compression device in market

F03M



Air Compression Therapy Pump

4/6/8 Chambers Available,
Massage A / Therapy B mode optional

Plug-in

F03001M

(Manual)



Easy Operation



Effective Lymphatic Edema

Specifications :
Device Size: 24cm*20cm*11cm
Mode Massage: A or Therapy B
Pressure Range: 30~250 mmHg
Treatment Time: 10~30 mins/always on
Power Consumption: 65VA

B Therapy Mode Sequential compression

F03001BM

(Manual)



Convenient One-Touch Sleeve



Use Various Cuffs According to User's Area

Specifications :
Device Size: 24cm*20cm*12cm
Mode Massage: A or Therapy B
Pressure Range: 30~250 mmHg
Treatment Time: 10~30 mins/always on
Power Consumption: 65VA

A Massage Mode

F03006

(Manual)



Convenient Handle Design

Home use
Specifications :
Device Size: 25cm*22cm*11cm
Mode Massage: A or Therapy B
Pressure Range: 30~250 mmHg
Treatment Time: 10~30 mins/always on
Power Consumption: 35VA

F03001

(Medical CE FDA510K)



Specifications :
Device Size: 24cm*20cm*11cm
Massage Mode: 6 modes
Pressure Range: 30 ~ 250 mmHg
Treatment Time: 10~ 60 mins
Power Consumption: 45VA
Each Chamber can be Partially Skipped
4,6,8 Chambers Available

F03008

(Medical FDA510K)



Specifications :
Device Size: 25cm*22cm*11cm
Massage Mode: 3 modes
Pressure Range: 30~250 mmHg
Treatment Time: 10~60 mins
Power Consumption: 65VA
Each Chamber can be Partially Skipped
4,6 Chambers Available



Use Durable Domestic Air Motors



6 Modes (Circulation & Massage)



Use Various Cuffs According to User's Area



A Choice of Sleeve



Set Pressure and Time



Skip Function



Effective Lymphatic Edema



Color Option

F03001B



Specifications :
Device Size: 24cm*20cm*12cm
Massage Mode: 6 modes
Pressure Range: 30 ~ 250 mmHg
Treatment Time: 10~ 60 mins
Power Consumption: 65VA
Each Chamber can be Partially Skipped
4,6,8 Chambers Available

F03001C



Specifications :
Device Size: 26.5cm*15.5cm*31cm
Massage Mode: 6 modes
Pressure Range: 30 ~ 250 mmHg
Treatment Time: 10~ 60 mins
Power Consumption: 65VA
Each Chamber can be Partially Skipped

With Li-ion Battery

F03002

(Medical CE FDA510K)



Specifications :

Device Size: 25cm*12cm*9cm
Massage Mode: 5 modes
Pressure Range: 30~250 mmHg
Treatment Time: 10~60 mins
Power Consumption: 65VA
Input: DC 15V/2600mA
Each Chamber can be Partially Skipped

F03002B



Specifications :

Device Size: 23cm*12.5cm*10cm
Massage Mode: 5 modes
Pressure Range: 30 ~ 250 mmHg
Treatment Time: 10~ 60 mins
Input: DC 15V/2600mA/5200mA
Each Chamber can be Partially Skipped

F03002C



Specifications :

Device Size: 27cm*12.5cm*10cm
Massage Mode: 5 modes
Pressure Range: 30 ~ 250 mmHg
Treatment Time: 10~ 60 mins
Input: DC 15V/2600mA/5200mA
Each Chamber can be Partially Skipped

F03002D



Specifications :

Device Size: 26.5cm*12.5cm*9cm
Massage Mode: 5 modes
Pressure Range: 30 ~ 250 mmHg
Treatment Time: 10~ 60 mins
Input: DC 15V/2600mA/5200mA
Each Chamber can be Partially Skipped

CE 0598 FDA

NEW

Hand Rehabilitation Device Stroke Patients Finger Recovery Training



F03002H



F03002G

Model	F03002H/F03002G
Device Size	25cm*12cm*9cm/14.5*12.5*5.5cm
Built-in Battery Volume	2600mAH (7.2V/12V)
Optional Programs	4 modes: A, B, C, D
Pressure Range	100mmHg~180mmHg
5 main functions	Single-finger, multi-finger training function, passive function, mirror function, air wave massage function.

Main Features

- Portable Compact Host Device, Low Mute.
- Multiple Functions: Single-Finger, Multi-Finger Training Function, Passive Function, Mirror Function, Air Wave Massage Function.
- Touch Keypad, Adjustable Treating Time with LED Display, and Available to Toggle Mode.
- The Recovery Glove for Finger Extension and Flexion is a Glove That Imitates the Hand Structure and Fits the Shape of the Hand, Based on The Pneumatic Working Principle, It Helps Hands Grip and Release Exercise.
- Wireless Data Glove Controls The Extension and Flexion of The Other Glove.
- Added Air Wave Massage Function, Built-in Pressure Sensor, Air Pressure Range is 100mmHg-180mmHg.
- Built 2600mA Lithium Battery Allows Longer Working Time, with a maximum use of 4 hours.

