

FOFO



www.fofomedical.com

Foshan Hongfeng Co.,Ltd

Address: No.4 Leqiang Road, Leping Sanshui, Foshan, 528100, GuangDong, China
 Tel: 0757-88791710 ext 888 / 0592-2097262
 Email: fofo@fofomedical.com
 Mob: +86 13055506651
 www.fofomedical.com

FOFO OVER **20** YEARS EXPERIENCE


IN HIGH-FREQUENCY PRODUCTS



COMPANY INTRODUCTION

Foshan Hongfeng Co., Ltd mainly produces medical rehabilitation products, specializing in the design and manufacture of medical air mattresses, water mattresses, air compression therapy machines, compressor nebulizers, suction machines, oxygen concentrators, blood pressure monitors, air cushions and other rehabilitation products. Foshan HongFeng has over 20 years of experience in making high-frequency products. Our mattress high frequency machines are 225KW, and are the most technologically advanced, energy-efficient and low-carbon machines in Mainland China. There are over 360 workers in the group. Our company has strong R & D capabilities, from mold design and manufacturing to software research and development, high-frequency product design and manufacturing are all completed internally, and quality control meets European and American requirements.

Our products have been approved by FDA, 510K, ISO and CE. The sales team's languages include English, Russian, Spanish, Portuguese, and Japanese. We have exported to the USA, Europe, Asia, the Middle East, and many other countries in the world. We have earned a great reputation in our industry amongst our customers and peers. We sincerely hope we can build a stable business relationship with you and welcome you to visit our factory. We guarantee high-quality products, competitive prices, prompt delivery schedules and reliable services.



CE FDA ISO13485





APPLICATION

Lymphedema treatment typically includes a few methods and compression devices is one of them. The pneumatic compression pump provides on-and-off pressure to keep fluid moving through lymph vessels and veins so it doesn't build up in arms, legs and other parts of your body. The device works by connecting the pump to a sleeve that wraps around affected areas. The pump drives a regular cycle of pressure that inflates and deflates the sleeve. Compression pumps may also reduce complications of lymphedema.

How can air compression therapy systems help?

- ✔ Improving blood and lymph circulation.
- ✔ Pushing the residual lymph towards the lymph nodes.
- ✔ Removal of lymphatic edema and exudates.
- ✔ Improving the nutrition and oxygenation of the skin and other tissues.



Avoid wound



Pain focus training

Air Compression Therapy System Cuffs

Sleeve	Description	Size
Leg	Leg Cuff M	91*65cm
	Leg Cuff L	100*70cm
	Leg Cuff XL	110*70cm Overlapping
	Leg Cuff XXXL	125*76cm Overlapping
Arm	Arm Cuff M	87*51cm
	Arm Cuff L	90*60cm Overlapping
Waist	Waist Cuff M	120*37cm
	Waist Cuff L	130*37cm Overlapping
Trousers	Short Pants (4,6,8 Chambers)	80*76cm
	Long Pants (6,8 Chambers)	142*85cm
Vest/cloth	Vest (6,8 Chambers)	78*72cm
Jacket	Double Arm/Shoulder	133*37 (pair)/cm



ABE DCF 6 Modes



A One by one, single part massage, after all chambers are treated, one cycle complete


B One after another, when all chambers are inflated together, one cycle complete

C Four times inflatable massage from bottom to top

D All together to inflate, and then deflate together, four times as one cycle

E Combined B and C mode


F Combined A and B mode




Intermittent Pressure Compression Device (IPC)

4/6/8 Chamber Medical CE FDA510K

- Powerful Medical-Grade Air Pump Provides 250mmHg Pressure
- Upgrad A-F 6 Modes
- Chambers Optional Function Skips Parts which you don't want to Squeeze
- Sleeves with Overlapping Design Provide more Comfort
- Extension Straps for Wider Sleeves
- Easy-Split Designs for 2-4 Sleeves Parallel use



Remote control distance <5 meters



FO-3001



Specifications:

Pump	Sleeves	
Item: FO3001, FO3008	Leg(4/6/8 Chambers)	M 91*65cm / L 100*70cm
Power Supply: 110V/220V/50-60Hz		XL 110*70cm overlapping / XXXL 125*76cm overlapping
Air Output: 15-20 liter/min.	Arm(4/6/8 Chambers)	M 87*51cm / L 90*60cm overlapping
Pressure Range:30mmHg - 250 mmHg	Waist(4/6/8 Chambers)	M 120*37cm / L 130*37cm overlapping
Time: 10-60 min.	Long Pants (6/8 Chambers)	142*85cm / 142*76cm overlapping
Power Consumption: 65VA	Short Pants (4/6/8 Chambers)	82*77cm
Size: 24cm(L) x 20cm(W) x 11cm(H)	Jacket (6/8 Chambers)	102*73cm
Weight: 2.45kg	Vest(6/8 Chambers)	78*72cm

Sequential Compression Device (SCD)

4/6/8 Chambers

- Medical CE, FDA510K ClassII
- Li-ion Battery 2600mA
- Powerful Medical-Grade Air Pump Provides 250mmHg Pressure
- Upgrade A-E 5 Modes Instead of the old Two-Modes Massage and Therapy

FO3002



Specifications:

Pump	Sleeves	
Item: FO3002	Leg(4/6/8 Chambers)	M 91*65cm / L 100*70cm
Power Supply: 100-220V		XL 110*70cm overlapping / XXXL 125*76cm overlapping
Air Output: 15-20 liter/min.	Arm(4/6/8 Chambers)	M 87*51cm / L 90*60cm overlapping
Pressure Range:30mmHg - 250 mmHg	Waist(4/6/8 Chambers)	M 120*37cm / L 130*37cm overlapping
Time: 10-60 min.	Long Pants (6/8 Chambers)	142*85cm / 142*76cm overlapping
Power Consumption: 65VA	Short Pants (4/6/8 Chambers)	82*77cm
Size: 25cm(L) x 12cm(W) x 9cm(H)	Jacket (6/8 Chambers)	102*73cm
Weight: 1-1.6 kg	Vest(6/8 Chambers)	78*72cm

NEW


4.3" Screen Air Compression Customer Therapy System

WiFi and Bluetooth Connectable/ App Control

SPECIFICATIONS

- Power Consumption: 30W/40VA
- Device Size: 25cm*12cm*9cm
- Massage Mode: 6 Modes
- Pressure Range: 30~250 mmHg
- Treatment Time: 10~60 mins
- Input: DC 24V/2600mA
- Each Chamber can be Partially Skipped

FO3002-1



Premium

LCD touch screen Air Compression Therapy System with Bluetooth



App Optional



FO3002-E

- 7-inch large screen
- APP synchronization control
- Pre set cold therapy mode
- 3-4 hours of battery life
- USB battery external discharge
- Habitual memory

WiFi and Bluetooth Connectable/ App Control

Built-in Video Instructions

Full-Body Garments


SPECIFICATIONS

- Multi-Lingual Device
- Demo Videos
- Power Consumption: 30W/40VA
- Device Size: 25cm*12cm*9cm
- Massage Mode: 6 modes
- Pressure Range: 30~270 mmHg
- Treatment Time: 10~60 mins
- Input: DC 15V/5200mA
- Each Chamber can be Partially Skipped
- Can use as USB Mobile Phone Charger

Massage mattress

Automatic Back Stretching Air pad, Lumbar Spine Relaxer, Lumbar Spine Stretcher, Waist Tensioner, Loose top, Lumbar Pressure, Back Exercises

FO3M



Heating/Cold+ Compression

Heating/Cryo/Pressure therapy

Pressurization and Heating Simultaneously Accelerate Blood Flow



SPECIFICATIONS

- Size: 23.5 x24x10.5cm
- Machine Voltage: AC 220V
- Power: 30W/40VA
- Weight: 2.2KG
- Heating Module voltage: Dc 15V
- Feeling Temperature: approx. 40°C



High Quality Air Compressor



Heated Leg Straps



Targeted Therapy for Foot, Calf, Knee or Thigh



Effective for Lymphatic Edema

Cold/Pressure Therapy

SPECIFICATIONS

- Model: FO5001
- Name: Cryo+Compression Therapy
- Dimension: 290*270*285mm
- Weight: 3.1KG
- Water Flow: 0.85L/MIN
- Volume: 6.5L
- Temperature: 5°C
- Full Leg Cooling and Compression



Air Compression Therapy Massager

4/6/8 Chamber Available,
Massage A / Therapy B mode optional

FO3001

(Manual)



Easy Operation



Effective Lymphatic Edema

Specifications :

- Device Size: 24cm*20cm*11cm
- Mode Massage: A or Therapy B
- Pressure Range: 30~250 mmHg
- Treatment Time: 10~60 mins/always on
- Power Consumption: 65VA

B Therapy Mode Sequential compression



FO3001B

(Manual)



Convenient One-Touch Sleeve



Use Various Cuffs According to User's Area

Specifications :

- Device Size: 24cm*20cm*12cm
- Mode Massage: A or Therapy B
- Pressure Range: 30~250 mmHg
- Treatment Time: 10~60 mins/always on
- Power Consumption: 65VA

A Massage Mode



FO3008

(Manual)



Convenient Handle Design

Home use

Specifications :

- Device Size: 25cm*22cm*11cm
- Mode Massage: A or Therapy B
- Pressure Range: 30~250 mmHg
- Treatment Time: 10~60 mins/always on
- Power Consumption: 65VA

Plug-in

FO-3001

(Medical CE FDA510K)



Specifications :

Device Size: 24cm*20cm*11cm
 Massage Mode: 6 modes
 Pressure Range: 30 ~ 250 mmHg
 Treatment Time: 10~ 60 mins
 Power Consumption: 65VA
 Each Chamber can be Partially Skipped
 4,6,8 Chambers Available

FO-3008

(Medical FDA510K)



Specifications :

Device Size: 25cm*22cm*11cm
 Massage Mode: 3 modes
 Pressure Range: 30~250 mmHg
 Treatment Time: 10~60 mins
 Power Consumption: 65VA
 Each Chamber can be Partially Skipped
 4,6 Chambers Available



FO-3001B

(Massage)



Specifications :

Device Size: 24cm*20cm*12cm
 Massage Mode: 6 modes
 Pressure Range: 30 ~ 250 mmHg
 Treatment Time: 10~ 60 mins
 Power Consumption: 65VA
 Each Chamber can be Partially Skipped
 4,6,8 Chambers Available

FO-3001C

(Massage)



Specifications :

Device Size: 26.5cm*15.5cm*31cm
 Massage Mode: 6 modes
 Pressure Range: 30 ~ 250 mmHg
 Treatment Time: 10~ 60 mins
 Power Consumption: 65VA
 Each Chamber can be Partially Skipped

With Li-ion Battery

FO3002

(Medical CE FDA510K)



Specifications :

Device Size: 25cm*12cm*9cm
 Massage Mode: 5 modes
 Pressure Range: 30~250 mmHg
 Treatment Time: 10~60 mins
 Power Consumption: 65VA
 Input: DC 15V/2600mA
 Each Chamber can be Partially Skipped

FO3002-B

(Massage)



Specifications :

Device Size: 23cm*12.5cm*10cm
 Massage Mode: 5 modes
 Pressure Range: 30 ~ 250 mmHg
 Treatment Time: 10~ 60 mins
 Input: DC 15V/2600mA/5200mA
 Each Chamber can be Partially Skipped

FO3002-C

(Massage)



Specifications :

Device Size: 27cm*12.5cm*10cm
 Massage Mode: 5 modes
 Pressure Range: 30 ~ 250 mmHg
 Treatment Time: 10~ 60 mins
 Input: DC 15V/2600mA/5200mA
 Each Chamber can be Partially Skipped

FO3002-D

(Massage)



Specifications :

Device Size: 26.5cm*12.5cm*9cm
 Massage Mode: 5 modes
 Pressure Range: 30 ~ 250 mmHg
 Treatment Time: 10~ 60 mins
 Input: DC 15V/2600mA/5200mA
 Each Chamber can be Partially Skipped

NEW Air compression leg massager

Main Features

- Improving blood circulation
- Prevent muscle atrophy and muscle fibrosis
- Physical recovery and muscles after exercise
- Improving the nutrition and oxygenation of the skin and other tissue
- 4 modes : massage and therapy function
- Pressure adjustable up to 160mmHg
- Ultra-light easy carry



FO3003B

Model	FO3003B
Device Size	20*6*5CM
Chamber	6 chamber sleeves
Modes	F1 individual massage /F2 sequential compression /F3 fully compression /F4 combined
Pressure	60/80/100/120/140/160mmHg
Power Consumption	AC100-230V,DC12V
Power	14W
Noise Level	Below 55 dBA
Certificates	massage CE



FO3002 MINI



Wireless Calf Mini Massager:

- Model No: FO3002min with Battery
- Mode: F1, F2, F3
- Leg Strap Size: 600*350*0.25cm
- Main Device Size: 14*7*0.25cm
- Pressure: 60/80/99mmHg
- Heating Temperature: 45°C-55°C
- Input Voltage DC5V/2200 mA
- Carton Box: 380*340*380mm; 10pcs per carton; 5800pcs for 20ft, 13870 for 40HQ
- Certificates: ISO, FDA, Massage CE

NEW Hand Rehabilitation Device
Stroke Patients Finger Recovery Training



FO3002H

Model	FO3002H
Device Size	25cm*12cm*9cm
Built-in Battery Volume	2600mAH (4 hours use)
Optional Programs	4 modes: A, B, C, D
Pressure Range	100mmHg~180mmHg
5 main functions	Single-finger, multi-finger training function, passive function, mirror function, air wave massage function.



FO3002G

Main Features

- Portable Compact Host Device, Low Mute.
- Multiple Functions: Single-Finger, Multi-Finger Training Function, Passive Function, Mirror Function, Air Wave Massage Function.
- Touch Keypad, Adjustable Treating Time with LED Display, and Available to Toggle Mode.
- The Recovery Glove for Finger Extension and Flexion is a Glove That Imitates the Hand Structure and Fits the Shape of the Hand, Based on The Pneumatic Working Principle, It Helps Hands Grip and Release Exercise.
- Wireless Data Glove Controls The Extension and Flexion of The Other Glove.
- Added Air Wave Massage Function, Built-in Pressure Sensor, Air Pressure Range is 100mmHg-180mmHg.
- Built 2600mA Lithium Battery Allows Longer Working Time, and The Longest can be Used for 4 Hours.

